


Filename: 2020-06-26-01_Monogenea_for_AMEB.gDNA


Gel Images



Default image (Contrast 100%), Image is Scaled to Sample

Sample Info


Well	DIN	Conc. [ng/μl]	Sample Description	Alert	Observations
A1	-	177	Ladder	⚠	Caution! Expired ScreenTape device; Ladder
B1	3.3	9.67	B1_54-55_Capsalidae_sp	⚠	Caution! Expired ScreenTape device; Sample concentration outside recommended range
C1	7.1	5.03	B1_56-60_Leptocotyle_minor	⚠	Caution! Expired ScreenTape device; Sample concentration outside recommended range
D1	5.7	23.6	B1_61-63_Dicycotyle_coeliaca	⚠	Caution! Expired ScreenTape device
E1	5.4	20.3	B1_64-65_1_Dicycotyle_coeliaca	⚠	Caution! Expired ScreenTape device
F1	4.9	10.1	B1_64-65_2_Dicycotyle_coeliaca_2	⚠	Caution! Expired ScreenTape device

A1: Ladder**Sample Table**

Well	DIN	Conc. [ng/μl]	Sample Description	Alert	Observations
A1	-	177	Ladder		Caution! Expired ScreenTape device; Ladder

Peak Table


Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	% Integrated Area	From [bp]	To [bp]	Peak Comment	Observations
100	8.50	8.50	-	66	150		Lower Marker
250	13.2	-	7.86	173	311		
400	14.0	-	8.37	325	481		
600	13.9	-	8.31	508	718		
900	13.4	-	7.99	786	1059		
1200	13.1	-	7.83	1059	1361		
1500	13.2	-	7.90	1361	1770		
2000	13.2	-	7.85	1770	2254		
2500	13.1	-	7.83	2254	2777		
3000	13.1	-	7.83	2777	3451		
4000	12.9	-	7.67	3587	5084		
7000	12.2	-	7.30	5809	10511		
15000	12.5	-	7.45	10511	20347		
48500	9.44	-	5.63	20347	>60000		
-	-	-	-	-	-		Sample Well

B1: B1_54-55_Capsalidae_sp**Sample Table**

Well	DIN	Conc. [ng/μl]	Sample Description	Alert	Observations
B1	3.3	9.67	B1_54-55_Capsalidae_sp		Caution! Expired ScreenTape device; Sample concentration outside recommended range

Peak Table


Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	% Integrated Area	From [bp]	To [bp]	Peak Comment	Observations
100	8.50	8.50	-	63	152		Lower Marker
1681	0.524	-	30.69	1607	1908		
10024	1.10	-	64.59	6589	33726		
-	-	-	-	-	-		Sample Well

C1: B1_56-60_Leptocotyle_minor**Sample Table**

Well	DIN	Conc. [ng/μl]	Sample Description	Alert	Observations
C1	7.1	5.03	B1_56-60_Leptocotyle_minor		Caution! Expired ScreenTape device; Sample concentration outside recommended range

Peak Table


Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	% Integrated Area	From [bp]	To [bp]	Peak Comment	Observations
100	8.50	8.50	-	64	154		Lower Marker
16505	2.63	-	91.92	5607	>60000		
-	-	-	-	-	-		Sample Well

D1: B1_61-63_Dicycotyle_coeliaca**Sample Table**

Well	DIN	Conc. [ng/μl]	Sample Description	Alert	Observations
D1	5.7	23.6	B1_61-63_Dicycotyle_coeliaca		Caution! Expired ScreenTape device

Peak Table

Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	% Integrated Area	From [bp]	To [bp]	Peak Comment	Observations
100	8.50	8.50	-	63	158		Lower Marker
15638	8.83	-	93.39	4728	>60000		
-	-	-	-	-	-		Sample Well


E1: B1_64-65_1_Dicycotyle_coeliaca**Sample Table**

Well	DIN	Conc. [ng/μl]	Sample Description	Alert	Observations
E1	5.4	20.3	B1_64-65_1_Dicycotyle_coeliaca		Caution! Expired ScreenTape device

Peak Table

Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	% Integrated Area	From [bp]	To [bp]	Peak Comment	Observations
100	8.50	8.50	-	65	145		Lower Marker
12141	8.43	-	98.17	3529	56158		
-	-	-	-	-	-		Sample Well

F1: B1_64-65_2_Dicycotyle_coeliaca_2



Sample Table

Well	DIN	Conc. [ng/μl]	Sample Description	Alert	Observations
F1	4.9	10.1	B1_64-65_2_Dicycotyle_coeliaca_2		Caution! Expired ScreenTape device

Peak Table

Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	% Integrated Area	From [bp]	To [bp]	Peak Comment	Observations
100	8.50	8.50	-	64	154		Lower Marker
12585	3.22	-	97.99	4215	56092		
-	-	-	-	-	-		Sample Well