

**Run Info**

Experiment Name	Dcoel010720
Sample ID	Dcoel
Run ID	87acb637-f3b9-4f14-b587-768f6b7270f1
Flow Cell Id	FAO01594
Start Time	July 1, 09:51
Run Length	7h 10m

**Run Summary**

Reads Generated	593.68 K
Bases Generated	571.21 Mb
Estimated Bases	516.87 Mb


**Run Parameters**

Flow Cell Type	FLO-MIN106
Kit	SQK-LSK109
Basecalling	on
Specified Run Length	72 hours
Initial Bias Voltage	-180 mV
FAST5 Output	Enabled
FAST5 Output Options	zlib_compress,fastq,raw
FASTQ Output Options	compress
FAST5 Reads per File	4000
FASTQ Output	Enabled
FASTQ Reads per File	4000
Active Channel Selection	Enabled
Mux Scan Period	1 hour 30 minutes
Reserved Pores	0 %
Basecall Model	Fast basecalling


**Versions**


MinKNOW Core	3.6.5
Bream	4.3.16
Guppy	3.2.10

Cumulative Output Reads





Cumulative Output Bases






The chart displays the State Time Equivalent (%) on the y-axis (0% to 100%) against Time on the x-axis. The data is represented by stacked bars with three colors: red (bottom), blue (middle), and grey (top). The overall trend shows a decrease in the State Time Equivalent over time, starting at approximately 73% at 5m and ending at approximately 48% at 6h 30m. There are several sharp drops in the value, notably at 1h 55m, 2h 50m, 4h 40m, and 6h 30m.


Time	Red (%)	Blue (%)	Grey (%)
5m	73	17	10
1h	88	12	0
1h 55m	63	37	0
2h 50m	56	44	0
3h 45m	73	27	0
4h 40m	49	51	0
5h 35m	58	42	0
6h 30m	40	50	10




Mux Scan Grouped




Mux Scan Categorised




Temperature History




Bias Voltage History




Translocation Speed



QScore



Disk Write Performance





## Run Debug Messages

- Flow cell FAO01594 has 591 pores available for sequencing. Starting sequencing with 318 pores July 1, 16:02
- Performing Mux Scan July 1, 16:00
- Flow cell FAO01594 has 727 pores available for sequencing. Starting sequencing with 372 pores July 1, 14:30
- Performing Mux Scan July 1, 14:28
- Flow cell FAO01594 has 892 pores available for sequencing. Starting sequencing with 413 pores July 1, 12:58
- Performing Mux Scan July 1, 12:57
- Flow cell FAO01594 has 1063 pores available for sequencing. Starting sequencing with 459 pores July 1, 11:27
- Performing Mux Scan July 1, 11:25
- Flow cell FAO01594 has 1263 pores available for sequencing. Starting sequencing with 497 pores July 1, 09:55
- Performing Mux Scan July 1, 09:53
- Starting sequencing procedure July 1, 09:53
- Waiting up to 300 seconds for temperature to stabilise at 34.0°C July 1, 09:51
- Disk / has 410 GB space remaining July 1, 09:51